

Corning[®] Supersomes[™] Ultra HEK293-derived Human Aldehyde Oxidase (AOX1) Cytosol

Catalog Number.....456801
Lot Number.....1067002

Storage Conditions.. -80°C
Date Released.....2021 September
Expiration Date.....2026 March

Product Description: Human AOX1 cytosol is prepared from Human Embryonic Kidney (HEK) 293 cells overexpressing human Aldehyde Oxidase (AOX1) (Gene Accession No. NM_001159). Aldehyde oxidase is a cytosolic enzyme that plays an important role in xenobiotic biotransformation. Primarily expressed in liver, AOX1 has broad substrate specificity, catalyzing the oxidation of aldehydes, aromatic azaheterocycles and nitrogen-carbon double bond containing compounds.

Negative Control Cytosol (Cat. No. 456800) prepared from HEK293 cells transfected with empty vector can be used as a control for this product.

Package Contents.....2.5 mg protein in 0.5 mL
Protein Content.....5.0 mg/mL in 0.1M potassium phosphate buffer, pH7.4

QUALITY CONTROL:

| Specification | Criteria | Result |
|-----------------|---|--------|
| Enzymatic Assay | Ultra Supersomes: Human AOX1 | PASS |
| | Substrate: 250 µM zaleplon | |
| | Incubation time: 60 min | |
| | Activity (pmol/mg/min) >=2 pmol/mg/min. | |

METHOD (zaleplon assay): The reaction is initiated by adding 10 µL of 5 mg/mL AOX1 cytosol into 90 µL 100 mM potassium phosphate (pH 7.4, 37°C) containing 278 µM zaleplon substrate. The reaction mixture is incubated for 60 min at 37°C. The reaction is terminated by adding 100 µL pre-chilled acetonitrile containing 0.4 µM labetalol as internal standard. The metabolite formation is quantified using LC-MS analysis.

STORAGE and HANDLING

- Thaw on ice and keep on ice until use.
- The activity can be compromised upon freezing and thawing. The activity remains 80% or above if the sample is thawed less than 4 times. After initial thaw, it is highly recommended to aliquot sample into single use vials to minimize freeze-thaw cycles.
- Triplicate samples per assay conditions are recommended.
- The concentration of test compound and incubation time usually depends upon the specific characteristics of the drug.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks


© 2017 Corning Incorporated

Discovery Labware, Inc.
6 Henshaw Street
Woburn, MA 01801
Tel: 1.978.442.2200 (U.S.)
CLSTechServ@Corning.com

CORNING

www.corning.com/lifesciences

Quality Assurance:

 _____

Date:

September 15, 2021 _____

SAFETY RECOMMENDATIONS:

Safety assessment indicates this product is not hazardous, therefore no SDS (Safety Data Sheet) is provided. Handle in accordance with good industrial hygiene and laboratory safety practices.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks

© 2017 Corning Incorporated