

## Corning® Supersomes™ Ultra HEK293-derived Human Aldehyde Oxidase (AOX1) Cytosol

Catalog Number.....456801  
Lot Number.....9046002

Storage Conditions.. -80°C  
Date Released.....2019 April  
Expiration Date.....2024 February

**Product Description:** Human AOX1 cytosol is prepared from Human Embryonic Kidney (HEK) 293 cells overexpressing human Aldehyde Oxidase (AOX1) (Gene Accession No. NM\_001159). Aldehyde oxidase is a cytosolic enzyme that plays an important role in xenobiotic biotransformation. Primarily expressed in liver, AOX1 has broad substrate specificity, catalyzing the oxidation of aldehydes, aromatic azaheterocycles and nitrogen-carbon double bond containing compounds.

Negative Control Cytosol (Cat. No. 456800) prepared from HEK293 cells transfected with empty vector can be used as a control for this product.

Package Contents.....2.5 mg protein in 0.5 mL  
Protein Content.....5.0 mg/mL in 0.1M potassium phosphate buffer, pH7.4

### QUALITY CONTROL:

Specification	Criteria	Result
Enzymatic Assay	Ultra Supersomes: Human AOX1	PASS
	Substrate: 250 µM zaleplon	
	Incubation time: 60 min	
	Activity (pmol/mg/min) >=2 pmol/mg/min.	

**METHOD (zaleplon assay):** The reaction is initiated by adding 10 µL of 5 mg/mL AOX1 cytosol into 90 µL 100 mM potassium phosphate (pH 7.4, 37°C) containing 278 µM zaleplon substrate. The reaction mixture is incubated for 60 min at 37°C. The reaction is terminated by adding 100 µL pre-chilled acetonitrile containing 0.4 µM labetalol as internal standard. The metabolite formation is quantified using LC-MS analysis.

### STORAGE and HANDLING

- Thaw on ice and keep on ice until use.
- The activity can be compromised upon freezing and thawing. The activity remains 80% or above if the sample is thawed less than 4 times. After initial thaw, it is highly recommended to aliquot sample into single use vials to minimize freeze-thaw cycles.
- Triplicate samples per assay conditions are recommended.
- The concentration of test compound and incubation time usually depends upon the specific characteristics of the drug.

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks)


© 2017 Corning Incorporated

Discovery Labware, Inc.  
6 Henshaw Street  
Woburn, MA 01801  
Tel: 1.978.442.2200 (U.S.)  
CLSTechServ@Corning.com

CORNING

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

Quality Assurance:

  
\_\_\_\_\_

Date:

28 April, 2019  
\_\_\_\_\_

**SAFETY RECOMMENDATIONS:**

Safety assessment indicates this product is not hazardous, therefore no SDS (Safety Data Sheet) is provided. Handle in accordance with good industrial hygiene and laboratory safety practices.

***For Research Use Only. Not for use in diagnostic or therapeutic procedures.***

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks)

© 2017 Corning Incorporated